



ATR KNOWS HOW TO HELP YOU GROW

**Large Capacity Incubator Shaker**  
**20-400 rpm Shaking Speed**  
**Available Options:**  
**CO<sub>2</sub> Control**  
**Refrigeration**  
**Humidity**  
**Lighting**  
**Slide-out Shelf**

## MULTITRON II



**INFORS KNOW-HOW PUTS YOU AHEAD**



Appropriate Technical Resources, Inc.  
9157 Whiskey Bottom Road, Laurel, MD 20723  
800-827-5931 [www.atrbiotech.com](http://www.atrbiotech.com)

The Multitron II is a large capacity incubator shaker. Multitron II may be configured as a single unit, or stacked as doubles or triples. The units allow for the optimum configuration in a limited space. The ingenious design of the INFORS shakers features a front-loading, pullout platform, stainless interior, a panoramic front window and inspection light. The patented drive provides uniform motion and quiet operation without vibration. The sealed cabinet ensures fast and accurate heating. Refrigeration, humidity and CO<sub>2</sub> control are available as optional controls.

## Operation

Each unit has a control panel with a splashproof keypad and LED display. A comfortable working height on each unit allows easy access to setpoints, on/off control and program cycles (day/night simulations). The display alternates between the current values of temperature and speed while the unit is operating.

Microprocessor technology and PID control guarantee the precise maintenance of setpoints. The machine will resume operation each time the door is closed or following a power failure. Whenever the door is opened, a brake is applied to the shaker drive and the circulation fans are switched off. The platform is released when the door is fully opened.

All input measurement signals are isolated and can be logged or used for external control via a RS232 serial link to a computer. Alarm modes are also available for external monitoring.

## Construction

The new Multitron II incubators each incorporate a fully independent control and power system. Each deck is fully redundant in terms of control of speed, temperature, lights, humidification, CO<sub>2</sub> and RS232 output. This new feature allows existing cabinets to be easily upgraded to a stacked system that saves valuable floor space. A single cabinet can be bench mounted or supplied with a fixed base support that allows for future expansion. All Multitron II Stacked systems come with the base support and stacking system.

The fully insulated incubation chamber coupled with the use of efficient blowers ensures accurate temperature control and minimal temperature variations over the surface of the platform.

Cooling can be fitted on any unit in a stacked system. Top or side mounted conventional refrigeration or rear mounted thermoelectric (Peltier) cooling can be specified depending on space limitations, ambient temperature and end temperatures requirements.

Lighting, CO<sub>2</sub>, and humidity, as options, can be installed in any cabinet location. Lighting requires cooling for stable temperature control.

The Multitron II can have a 25 mm or 50 mm diameter orbit (or throw) to optimize the oxygen transfer rate for different sized flasks.

Speed control ranges between 20 and 400 rpm. Cabinet location, throw, and the size of flasks determine maximum speed.

The cabinet base is formed to be a water tight tray and comes equipped with a drain fitting for complete containment of spills. The drive motor is located outside of and above the environmental cabinet for reliability and safety.

A work light is provided on the cabinet. Activated by pressing any key on the front panel, this allows you to see the contents of the shaker without opening the door.



With a modest cabinet footprint the MULTITRON II platform (850 x 470mm) accepts large numbers of Erlenmeyer flasks. The Universal and dedicated platforms are designed with threaded holes for mounting of INFORS clamps and test tube racks. The slide-out shelf option allows the Multitron II to double as a static incubator. The shelf is easy to load and is also removable. When installed, the slide-out shelf limits the flask height to 12".



As an alternative to the traditional clamp style platform INFORS offers an adhesive material platform of "Green Sticky Stuff" for maximum flexibility of different sized flasks and carriers. "Green Sticky Stuff", with its proprietary formula, is remarkable in its ability to hold flat bottom vessels. It can be rejuvenated by simple washing with water and liquid dish detergent. Customers can request combination platforms of "Green Sticky Stuff", flask clamps and tube racks.

## External Dimensions and Weight

Single Unit	21" (H) x 42" (W) x 34" (D)	Wt. 195 lbs.
Double stack	46" (H) x 42" (W) x 34" (D)	Wt. 390 lbs.
Triple stack	67" (H) x 42" (W) x 34" (D)	Wt. 585 lbs.
Triple w/ top cooling	78" (H) x 42" (W) x 34" (D)	Wt. 695 lbs.
Units w/ side cooling	as above except 53" (W)	add 40 lbs / cabinet

## Specifications

<b>Drive:</b>	Brushless AC induction motor with drive belt
<b>Shaking Speed:</b>	20 to 400 rpm
<b>Orbital Throw:</b>	25 or 50 mm diameter
<b>Accuracy:</b>	± 1% at maximum speed
<b>Temperature Range :</b>	5°C above ambient to 65°C without cooling 15°C below ambient to 65°C with cooling
<b>Accuracy:</b>	± 0.2 °C (Pt 100 sensor)
<b>Temperature Uniformity:</b>	± 0.5 °C over working surface
<b>Air Circulation:</b>	210 CFM (360 m <sup>3</sup> /h)
<b>Data Interface:</b>	RS232 bi-directional
<b>Power Supply:</b>	115V @ 60 Hz standard cabinet 9 Amps w/ cooling & lighting 13 Amps

Models		Individual Options	
Master Unit, 25 mm	AJ118	<b>Cooling System</b>	
Master Unit, 50 mm	AJ119	Cooling in side box	ZK120
<b>Base</b>		Peltier cooling	ZKPC
Trolley	ZM101	Casing w/cooling	ZM110
Support 6"	ZM106	<b>Lighting/cabinet</b>	
Support 9"	ZM109	Daylight	ZL110
<b>Cover Frame</b>	ZM110	Photosynthesis	ZL111
<b>Fixing Bars</b>		Algae	ZL116
Anchor Bar double	ZM111	<b>Slide-out shelf</b>	SS110
Anchor Bar triple	ZM110	<b>Humidity Control</b>	HC110
		<b>CO<sub>2</sub> Control</b>	CO110

Fixed Trays and Clamps		Additional Accessories	
<b>Universal with threaded holes</b>	<b>ATM0</b>	<b>Universal with sticky stuff</b>	<b>ATM00</b>
For fitting different clamps, test tubes racks, etc.		For flasks, bottles or any container with a wide, flat base.	
<b>Fixed Trays</b>	<b>ATM...</b>	<b>Universal Tray for Microtiter plates</b>	<b>ATM30</b>
Fitted with clamps for one size of Erlenmeyer flasks.			
<b>Fixed Trays with Stainless Steel Clamps</b>			
<b>ATM50</b>	50 ml	131 clamps	
<b>ATM100</b>	100 ml	87 clamps	
<b>ATM250</b>	250 ml	48 clamps	
<b>ATM500</b>	500 ml	31 clamps	
<b>ATM1000</b>	1000 ml	18 clamps	
<b>ATM2000</b>	2000 ml	13 clamps	
<b>ATM4000</b>	4000 ml	8 clamps	
<b>Additional Clamps</b>			
<b>ATK25</b>		25 ml clamp	
<b>ATK50</b>		50 ml clamp	
<b>ATK100</b>		100 ml clamp	
<b>ATK250</b>		250 ml clamp	
<b>ATK500</b>		500 ml clamp	
<b>ATK1000</b>		1000 ml clamp	
<b>ATK2000</b>		2000 ml clamp	
<b>ATK3000</b>		3000 ml clamp	
<b>ATK4000</b>		4000 ml clamp	
<b>ATK5000</b>		5000 ml clamp	
<b>Tube Racks</b>			
<b>ATH001</b>		1 roll Sticky Stuff 5 m	
<b>ATH008</b>		8 rolls Sticky Stuff 5 m	
<b>ATKF2</b>		1.8L Fernbach	
<b>ATKF3</b>		2.8L Fernbach	
<b>ATR 08</b>		Tube rack 144 x 8 mm	
<b>ATR 10</b>		Tube rack 100 x 10 mm	
<b>ATR 12</b>		Tube rack 72 x 12 mm	
<b>ATR 14</b>		Tube rack 64 x 14 mm	
<b>ATR 16</b>		Tube rack 42 x 16 mm	
<b>ATR 18</b>		Tube rack 39 x 18 mm	
<b>ATR 20</b>		Tube rack 36 x 20 mm	
<b>ATR 22</b>		Tube rack 33 x 22 mm	
<b>ATR 25</b>		Tube rack 20 x 25 mm	
<b>ATR 30</b>		Tube rack 15 x 30 mm	