

TRIPLE SANDWICH IMMUNOPRECIPITATION SYSTEM

Ag\* -mMab- rat anti-mouse Kappa chain MAb-Staph.aureus

GENERAL:

This is an indirect immunoprecipitation system for use with mouse MAb. Mouse immunoglobulin subclasses bind to protein A on the surface of Staph. aureus bacteria through their Fc portion. The IgG's vary in their capability to bind depending on the heavy chain subclass. For example: IgG2a and IgG2b bind better than some IgG1.

Mouse MAb contain K light chain  $\geq$  90%. A rat anti-mouse Kappa chain MAb (187-1) [Yelton, D. E., C. Desaynard and M. D. Scharff, 'Use of monoclonal anti-mouse immunoglobulin to detect mouse antibodies'; Hybridoma 1: 5-11, 1981] is a useful reagent for most mouse MAb independently of their Ig subclass. This reagent is a rat IgG2c and binds to protein A. Excellent results can be obtained with the combination of both 187-1 MAb and Staph. aureus.

METHOD:1) Preclearing of labeled cell lysates with Staph. aureus (SA).

Before immunoprecipitation, labeled cell lysates in a 1.5 ml Sarstedt conical tube should be mixed and incubated with Staph. aureus bacteria (200  $\mu$ l of a 10% suspension of SA per ml of cell lysates) to preclear the proteins that bind nonspecifically to SA. After shaking for 1 h at 4 C, the mixture is centrifuged (2,000 rpm x 10 min.) and the supernatant can be used for immunoprecipitation.

Before use, SA is washed 3 times with and resuspended in TSA pH 7.4, 1% TX-100 and 0.1% glycerol.

2) Immunoprecipitation. Precoat 1.5 ml Sarstedt tubes without caps with 1% Hb in TSA. Add in order wash buffer (25  $\mu$ l), 50-100  $\mu$ l of mouse MAb culture supernatant (the equivalent to 10  $\mu$ g)<sub>1,2,3</sub> and mix with 10-100  $\mu$ l labeled cell lysates containing 2-3 x 10<sup>6</sup> cpm of I-antigen. Vortex and allow to stand 2 h, at 4 C. Then, add 100  $\mu$ l/sample of rat anti-mouse Kappa chain 187-1 MAb culture supernatant ( $\approx$  7.4  $\mu$ g), vortex and allow to stand at 4 C for 0.5 h or overnight. After the second incubation 50  $\mu$ l of a 10% suspension of washed Staph. aureus per sample is added. Shake for 10 min. at 4 C. Washing and later steps are as described in in 'Immunoprecipitation for Analysis of Radiolabeled Cell Surface Molecules.'